

Counter electrode, Platinum *No. 04-00002*

Thank you for choosing a Sensolytics miniaturized counter electrode.

The Sensolytics counter electrode consists of a platinum coil imbedded in a stabilizing glass body and is ready to use.

Connection with the potentiostat is realized via the gold connection pin at the rear side of the counter electrode. Please avoid traction and tension to the connection pin to ensure a reliable electrical contact.

To clean the Platinum coil, rinse the electrode alternating with water and if required with ethanol (or acetone, isopropanol etc.).